


ORAL SYRINGES

CLASSIC




CUSTOMIZATION

- Marking on plunger
- Graduations (ml, mg, kg, etc.)
- Individual wrapping on request
- Accessories (syringe holder, syringe wiper, protective sheath)

BENEFITS

- Compliance with international regulations, including but not limited to FDA (listing in Drug Master File) and CE marking
- Suitable for pediatric medicines, suspensions, or drinkable solutions
- Wide capacity range
- Easy-to-use
- Oral syringe plunged into the bottle



A length
B length to collar

	CLASSIC SYRINGES	REFERENCE	MATERIALS Barrel / Plunger Resin	LENGTH	LENGTH TO COLLAR	MARKING
1	0.3 ml	CLA0.3	PE/PS	76 mm	67 mm	Plunger
2	0.6 ml	CLA0.6	PE/PS	112 mm	101 mm	Plunger
3	1.5 ml	CLA1.5	PE/PS	80 mm	70 mm	Plunger
4	C2 ml	CLAC2	PE/PS	103 mm	95 mm	Plunger
5	2 ml	CLA2	PE/PS	109 mm	98 mm	Plunger
6	TH5 ml	CLATH5	PE/PS	100 mm	90 mm	Plunger
7	5 ml	CLA5	PE/PS	121 mm	100 mm	Plunger
8	C5 ml	CLAC5	PE/PS	123 mm	112 mm	Plunger
9	6 ml	CLA6	PE/PS	136 mm	125 mm	Plunger
10	8 ml	CLA8	PE/PS	121 mm	104 mm	Plunger
11	10 ml	CLA10	PE/PS	155 mm	143 mm	Plunger




CUSTOMIZATION

- Marking on plunger or barrel
- Graduations (ml, mg, kg, etc.)
- Individual wrapping on request
- Accessories (syringe holder, PIBA)

BENEFITS

- Compliance with international regulations, including but not limited to FDA (listing in Drug Master File) and CE marking
- Suitable for pediatric medicines, suspensions, or drinkable solutions
- Strong legibility and easy administration
- Easy-to-use : dosage and administration, high legibility, low gliding force, long lasting
- Easily recyclable
- Plunged into the bottle or with a Press-In Bottle Adaptor (PIBA)
- Very hygienic when used with a Press-In Bottle Adaptor (PIBA)



A length
B length to collar

	ADVANCED SYRINGE	REFERENCE	MATERIALS Barrel / Plunger Resin	LENGTH	LENGTH TO COLLAR	MARKING
1	1 ml PF	ADV1PF	PP/PE	98.2 mm	86 mm	Barrel
2	3 ml PF	ADV3PF	PP/PE	88.85 mm	79.80 mm	Barrel/plunger
3	TH5 ml LS	ADVTH5LS	PP/PE	100 mm	91 mm	Barrel/plunger
4	TH5 ml PF	ADVTH5PF	PP/PE	108 mm	98 mm	Barrel/plunger
5	C10 ml PF	ADVC10PF	PP/PE	114 mm	104 mm	Barrel/plunger
6	C10 ml	ADVC10	PP/PE	114 mm	104 mm	Barrel/plunger




CUSTOMIZATION

- Marking on plunger or barrel Graduations (ml, mg, kg, etc.)
- Individual wrapping on request
- Accessories (syringe holder, PIBA)

BENEFITS

- Compliance with international regulations, including but not limited to FDA (listing in Drug Master File) and CE marking
- Suitable for challenging substances (alcohol, viscous solutions)
- Easy-to-use : dosage and administration, high legibility, low gliding force, long lasting
- Excellent seal for watertightness and dosing



A length
B length to collar

	EVOLUTION SYRINGE	REFERENCE	MATERIALS Barrel / Plunger Resin / Oring resin	LENGTH	LENGTH TO COLLAR	MARKING
1	1 ml PF JT	EVO1PF	PP/PE/silicone	98 mm	86 mm	Barrel
2	2 ml ORP	EVO2ORP	PP/PP/silicone	110 mm	97 mm	Barrel
3	3 ml PF JT	EVO3PF	PP/PP/silicone	90.2 mm	79.5 mm	Barrel/plunger
4	TH5 ml PF JT	EVOTH5PF	PP/PP/silicone	110 mm	98 mm	Barrel/plunger
5	10 ml PF JT	EVOC10PF	PP/PP/silicone	114 mm	104 mm	Barrel/plunger

ORAL SYRINGES

ACCESSORIES

Our accessories for oral syringes include syringe wipers, protective sheath, Press-In Bottle Adaptor (PIBA) and syringe holders.



ORAL SYRINGE		1 ORAL SYRINGE WIPER (PP RESIN)	2 PROTECTIVE SHEATH (PP RESIN)	CE marked 3 PIBA (PE RESIN) PP28 GLASS OR PLASTIC BOTTLE	4 PIBA + CAP CRC CAP PP28 PIBA GLASS BOTTLE TE CAP DIN18 PIBA GLASS BOTTLE	5 SYRINGE HOLDER (PP RESIN)
Classic	0.3 ml					
	0.6 ml					X
	1.5 ml	X	X			X
	C2 ml	X				X
	2 ml		X			X
	TH5 ml					X
	5 ml	X	X			X
	C5 ml					X
	6 ml	X	X			X
	8 ml					
	10 ml					
Advanced	1 ml PF			X	X	X
	3 ml PF			X	X	X
	TH5 ml LS					X
	TH5 ml PF			X	X	X
	C10 ml PF			X	X	X
	C10 ml					
Evolution	1 ml PF JT			X	X	X
	2 ml ORP					X
	3 ml PF JT			X	X	X
	TH5 ml PF JT			X	X	X
	10 ml PF JT			X	X	X