


ORAL SYRINGES CLASSIC




CUSTOMIZATION

- Marking on plunger
- Graduations (ml, mg, kg, etc.)
- Individual wrapping on request
- Accessories (syringe holder, syringe wiper, protective sheath)

BENEFITS

- Compliance with international regulations, including but not limited to FDA (listing in Drug Master File) and CE marking
- Suitable for pediatric medicines, suspensions, or drinkable solutions
- Wide capacity range
- Easy-to-use
- Oral syringe plunged into the bottle



A length
B length to collar

	CLASSIC SYRINGES	REFERENCE	MATERIALS Barrel / Plunger Resin	LENGTH	LENGTH TO COLLAR	MARKING
1	0.3 ml	CLA0.3	PE/PS	76 mm	67 mm	Plunger
2	0.6 ml	CLA0.6	PE/PS	112 mm	101 mm	Plunger
3	1.5 ml	CLA1.5	PE/PS	80 mm	70 mm	Plunger
4	C2 ml	CLAC2	PE/PS	103 mm	95 mm	Plunger
5	2 ml	CLA2	PE/PS	109 mm	98 mm	Plunger
6	TH5 ml	CLATH5	PE/PS	100 mm	90 mm	Plunger
7	5 ml	CLA5	PE/PS	121 mm	100 mm	Plunger
8	C5 ml	CLAC5	PE/PS	123 mm	112 mm	Plunger
9	6 ml	CLA6	PE/PS	136 mm	125 mm	Plunger
10	8 ml	CLA8	PE/PS	121 mm	104 mm	Plunger
11	10 ml	CLA10	PE/PS	155 mm	143 mm	Plunger

ORAL SYRINGES ADVANCED



CUSTOMIZATION

- Marking on plunger or barrel
- Graduations (ml, mg, kg, etc.)
- Individual wrapping on request
- Accessories (syringe holder, PIBA)

BENEFITS

- Compliance with international regulations, including but not limited to FDA (listing in Drug Master File) and CE marking
- Suitable for pediatric medicines, suspensions, or drinkable solutions
- Easy-to-use : dosage and administration, high legibility, low gliding force, long lasting
- Easily recyclable
- Plunged into the bottle or used with a Press-In Bottle Adaptor (PIBA)
- Very hygienic, low contamination, reduced access to the drug when used with a Press-In Bottle Adaptor (PIBA)



	ADVANCED SYRINGE	REFERENCE	MATERIALS Barrel / Plunger Resin	LENGTH	LENGTH TO COLLAR	MARKING
1	1 ml PF	ADV1PF	PP/PE	98.2 mm	86 mm	Barrel
2	3 ml PF	ADV3PF	PP/PE	88.85 mm	79.80 mm	Barrel/plunger
3	TH5 ml LS	ADVTH5LS	PP/PE	100 mm	91 mm	Barrel/plunger
4	TH5 ml PF	ADVTH5PF	PP/PE	108 mm	98 mm	Barrel/plunger
5	C10 ml PF	ADVC10PF	PP/PE	114 mm	104 mm	Barrel/plunger
6	C10 ml	ADVC10	PP/PE	114 mm	104 mm	Barrel/plunger

ORAL SYRINGES EVOLUTION




CUSTOMIZATION

- Marking on plunger or barrel
- Graduations (ml, mg, kg, etc.)
- Individual wrapping on request
- Accessories (syringe holder, PIBA)

BENEFITS

- Compliance with international regulations, including but not limited to FDA (listing in Drug Master File) and CE marking
- Easy-to-use : dosage and administration, high legibility, low gliding force, long lasting
- Excellent seal for watertightness and dosing
- Plunged into the bottle or used with a Press-In Bottle Adaptor (PIBA)
- Very hygienic, low contamination, reduced access to the drug when used with a Press-In Bottle Adaptor (PIBA)



A length
B length to collar

	EVOLUTION SYRINGE	REFERENCE	MATERIALS Barrel / Plunger Resin / Oring resin	LENGTH	LENGTH TO COLLAR	MARKING
1	1 ml PF JT	EVO1PF	PP/PE/silicone	98 mm	86 mm	Barrel
2	2 ml ORP	EVO2ORP	PP/PP/silicone	110 mm	97 mm	Barrel
3	3 ml PF JT	EVO3PF	PP/PP/silicone	90.2 mm	79.5 mm	Barrel/plunger
4	TH5 ml PF JT	EVOTH5PF	PP/PP/silicone	110 mm	98 mm	Barrel/plunger
5	10 ml PF JT	EVOC10PF	PP/PP/silicone	114 mm	104 mm	Barrel/plunger

ORAL SYRINGES ACCESSORIES



ORAL SYRINGE	1	2	3	4		5
	ORAL SYRINGE WIPER (PP RESIN)	PROTECTIVE SHEATH (PP RESIN)	PIBA (PE RESIN) PP28 GLASS OR PLASTIC BOTTLE	CAP WITH PIBA INTEGRATED CRC CAP PP28 PIBA GLASS BOTTLE	TE CAP DIN18 PIBA GLASS BOTTLE	SYRINGE HOLDER (PP RESIN)
Classic	0.3 ml					
	0.6 ml					X
	1.5 ml	X	X			X
	C2 ml	X				X
	2 ml		X			X
	TH5 ml					X
	5 ml	X	X			X
	C5 ml					X
	6 ml	X	X			X
	8 ml					
	10 ml					
Advanced	1 ml PF			X	X	X
	3 ml PF			X	X	X
	TH5 ml LS					X
	TH5 ml PF			X	X	X
	C10 ml PF			X	X	X
	C10 ml					
Evolution	1 ml PF JT			X	X	X
	2 ml ORP					X
	3 ml PF JT			X	X	X
	TH5 ml PF JT			X	X	X
	10 ml PF JT			X	X	X