



This innovative asthma **spacer** is specially made for dispensing inhaled treatment. Its new design is user friendly and permits a more comfortable dispensation for the child or for the care giver. And the most important, its efficiency in administering the treatment is higher.

BENEFITS

- Made in Italy and France
- Antistatic material
- Flexible and quick insertion of mostly all MDI
- Rotating handle: adjust the inclination of MDI to the child position and adapt to both right or left handed patient
- 2 valves to monitor the efficiency
- Fully dismantled for easy cleaning and drying
- Specific shape: increase the fine particles distribution, the patient breathes 10 times more fine particles
- CE marking (Class I medical device)