

# DROPPERS

## EAR DROPPERS



- 1 **Cap:** screw or clip-on
- 2 **Soft body** (BA or BX) or cannula. Body of BX model snapped on bottle neck size WI20
- 3 **Closure (tamper-evident or not):** attaches the dropper to the bottle

### CUSTOMIZATION

- Individual packaging in soft blister on request

### BENEFITS

- Compliance with international regulations, including but not limited to FDA and CE marking (according to the references)
- Administration of calibrated drops
- Soft and tapered tip to ease administration
- Closures available to complement the range

### EAR DROPPERS BA

REF.	BOTTLE NECK SIZE	DROP SIZE IN $\mu\text{L}$	BODY MATERIALS	CAP MATERIALS	RING MATERIALS
XBA42111200-DIVER	DIN18	50	PVC WHITE	PE WHITE	PP WHITE
XBA540412101	DIN18	55	EVA NATURAL		
XBA12112200-DIVER	PH18	50	PVC WHITE		
XBA12115200-DIVER	PH20	50	PVC WHITE		
XBA34044200-DIVER	PP28	55	EVA NATURAL		

### EAR DROPPERS BX

REF.	BOTTLE NECK SIZE	DROP SIZE IN $\mu\text{L}$	BODY MATERIALS	CAP MATERIALS	RING MATERIALS
XBX1201100-DIVER	WI20	40	PVC NATURAL	PE WHITE	NA
XBX4212110-DIVER	WI20	50	PVC NATURAL		
XBX4212100-DIVER	WI20	50	PVC WHITE		

### EAR DROPPERS - CANNULA

REF.	BOTTLE NECK SIZE	DROP SIZE IN $\mu\text{L}$	BODY MATERIALS	CAP MATERIALS	RING MATERIALS
32106NM	20-410	NA	PVC NATURAL	PE WHITE	PP WHITE
	20-415				
32107	24-410				
32108	24-415				

### DROPPER CANNULA REFILLS

CANNULA + CAP	REFERENCE	CANNULA RESIN	CANNULA $\emptyset$	CAP RESIN
	32402	PVC	8 mm	PE
	32401		18 mm	
	32400		20 mm	