

IQ-DOSE

INTUITIVE & QUICK DOSE

Easy to use, the innovative **IQ-Dose solution** is used to count and administer mini-tablets at the same time, without touching them. This new drug delivery format is ideal for adjusting the dose to the patient's need thanks to its high concentration of active ingredients. **IQ-Dose** is the perfect solution for clinical trials, pediatric medicine, Parkinson's disease, oncology, psychiatric treatments, etc. and offers complete hygiene and safety.

BENEFITS

- Select the exact quantity of mini-tablets prescribed at the same time
- Adjust the dose gradually depending on the response to treatment
- Hygienic and accurate: no need to count out mini-tablets one by one, reducing the risks of error and spillage
- Customizable: adaptable to different sizes and quantity of tablet, sizes of container and logo
- Compliance with international regulations, including but not limited to FDA and CE marking



Model presented
1 to 24 mini-tablets



X'TEMP[®]



X'temp is a system to reconstitute medicine or food supplements just before the administration. The instant mixture of active principles to be added at the last moment can be in liquid or powder form.

It is hermetically stored in the reservoir in the cap until it is tapped, when it falls into the liquid in the bottle to be mixed ready for use.

Ideal for food supplements or energy products, **X'temp** is also perfect for medicines.

It consists of a perforating cap with a tamper-evident ring and sealed, single-dose container.

BENEFITS

- Practical and easy to use: reconstitute medicines without opening the cap
- Fewer preservatives needed: because the two phases are separated until the last moment, they remain stable for longer ensuring the perfect medicine integrity
- Safe, reliable system: lid with tamper-evident ring, with the powder or liquid in a hermetically sealed reservoir

PRESS&TAKE

Intuitive and accurate, **Press&Take** delivers the right dose with a simple press.

Its reliably delivers irregularly shaped solid forms of medicine without breaking or getting them trapped in the mechanism.

A tablet is released on each press!

BENEFITS

- Easy-to-use
- Tolerates irregularities in shape/diameter of tablets

CUSTOMIZATION

- Adaptable to different sizes of pellets
- Customized marking and color



CAPSU'MATIC

Hygienic, safe and effective, **Capsu'Matic** dispenses all forms of pills or gel capsules, one-by-one, with a simple press. Only the capsules that are about to be taken quit the container, leaving the rest hygienically inside.

USE

- Food supplements, vitamins and medicines with high added value, where dosing or toxicity needs to be controlled
- Gel capsule sizes 00 to 3
- Two standard container sizes available: 30 and 60 gel capsules

BENEFITS

- Hygienic and ergonomic: only the medicine being taken leaves the container
- Secure opening/closing prevents dispersion of gel capsules
- Ecological: no need for outer wrapping of secondary packaging
- No need for laboratories to change existing manufacturing processes

CUSTOMIZATION

- Adaptable to different shapes of gel capsule, tablet, etc.
- Size, color, labelling
- Child safety and tamper-evident labelling
- Optional space for dessicant



This innovative asthma **spacer** is specially made for dispensing inhaled treatment. Its new design is user friendly and permits a more comfortable dispensation for the child or for the care giver. And the most important, its efficiency in administering the treatment is higher.



BENEFITS

- Made in Italy and France
- Antistatic material
- Flexible and quick insertion of mostly all MDI
- Rotating handle: adjust the inclination of MDI to the child position and adapt to both right or left handed patient
- 2 valves to monitor the efficiency
- Fully dismantled for easy cleaning and drying
- Specific shape: increase the fine particles distribution, the patient breathes 10 times more fine particles
- CE marking (Class I medical device)

DISPENSING PACIFIER

With the Dispensing Pacifier, administering liquid is no longer a fight.

The dispensing pacifier is a comfortable way to give young patients any prescriptions they need, with no added stress for parents. Perfectly adapted for pediatric use, this device makes it a lot easier to administer liquid supplements and treatments.

The pacifier consists of five pieces: a nipple with a reservoir and a protective cap, a bracket and a link.

USE

1. Unscrew the reservoir and fill it with the liquid according to the prescribed dose
2. Screw the reservoir on the nipple bracket
3. Remove the protective cover and put the pacifier into the baby's mouth

