


ORAL SYRINGES CLASSIC




CUSTOMIZATION

- Marking on plunger
- Graduations (ml, mg, kg, etc.)
- Individual wrapping on request
- Accessories (syringe holder, syringe wiper, protective sheath)

BENEFITS

- Compliance with international regulations, including but not limited to FDA (listing in Drug Master File) and CE marking
- Suitable for pediatric medicines, suspensions, or drinkable solutions
- Wide capacity range
- Easy-to-use
- Oral syringe plunged into the bottle



A length
B length to collar

	CLASSIC SYRINGES	REFERENCE	MATERIALS Barrel / Plunger Resin	LENGTH	LENGTH TO COLLAR	MARKING
1	0.3 ml	CLA0.3	PE/PS	76 mm	66 mm	Plunger
2	0.6 ml	CLA0.6	PE/PS	112 mm	101 mm	Plunger
3	1.5 ml	CLA1.5	PE/PS	78 mm	69 mm	Plunger
4	2 ml	CLAC2	PE/PS	102 mm	94 mm	Plunger
5	2 ml	CLA2	PE/PS	109 mm	97 mm	Plunger
6	5 ml	CLATH5	PE/PS	100 mm	90 mm	Plunger
7	5 ml	CLA5	PE/PS	121 mm	100 mm	Plunger
8	5 ml	CLAC5	PE/PS	123 mm	112 mm	Plunger
9	6 ml	CLA6	PE/PS	136 mm	125 mm	Plunger
10	8 ml	CLA8	PE/PS	121 mm	103 mm	Plunger
11	10 ml	CLA10	PE/PS	154 mm	143 mm	Plunger

ORAL SYRINGES ADVANCED




CUSTOMIZATION

- Marking on plunger or barrel
- Graduations (ml, mg, kg, etc.)
- Individual wrapping on request
- Accessories (syringe holder, PIBA)

BENEFITS

- Compliance with international regulations, including but not limited to FDA (listing in Drug Master File) and CE marking
- Suitable for pediatric medicines, suspensions, or drinkable solutions
- Easy-to-use : dosage and administration, high legibility, low gliding force, long lasting
- Easily recyclable
- Plunged into the bottle or used with a Press-In Bottle Adaptor (PIBA)
- Very hygienic, low contamination, reduced access to the drug when used with a Press-In Bottle Adaptor (PIBA)




A length
B length to collar

	ADVANCED SYRINGE	REFERENCE	MATERIALS Barrel / Plunger Resin	LENGTH	LENGTH TO COLLAR	MARKING
1	1 ml PF	ADV1PF	PP/PE	98 mm	86 mm	Barrel
2	3 ml PF	ADV3PF	PP/PE	89 mm	80 mm	Barrel/plunger
3	TH5 ml LS	ADVTH5LS	PP/PE	100 mm	91 mm	Barrel/plunger
4	TH5 ml PF	ADVTH5PF	PP/PE	108 mm	98 mm	Barrel/plunger
5	C10 ml PF	ADVC10PF	PP/PE	114 mm	104 mm	Barrel/plunger
6	C10 ml	ADVC10	PP/PE	114 mm	104 mm	Barrel/plunger

ORAL SYRINGES

EVOLUTION




CUSTOMIZATION

- Marking on plunger or barrel
- Graduations (ml, mg, kg, etc.)
- Individual wrapping on request
- Accessories (syringe holder, PIBA)

BENEFITS

- Compliance with international regulations, including but not limited to FDA (listing in Drug Master File) and CE marking
- Easy-to-use : dosage and administration, high legibility, low gliding force, long lasting
- Excellent seal for watertightness and dosing
- Plunged into the bottle or used with a Press-In Bottle Adaptor (PIBA)
- Very hygienic, low contamination, reduced access to the drug when used with a Press-In Bottle Adaptor (PIBA)



	EVOLUTION SYRINGE	REFERENCE	MATERIALS Barrel / Plunger Resin / Oring resin	LENGTH	LENGTH TO COLLAR	MARKING
1	1 ml PF JT	EVO1PF	PP/PE/silicone	98 mm	86 mm	Barrel
2	2 ml ORP	EVO2ORP	PP/PP/silicone	110 mm	97 mm	Barrel
3	3 ml PF JT	EVO3PF	PP/PP/silicone	90 mm	80 mm	Barrel/plunger
4	TH5 ml PF JT	EVOTH5PF	PP/PP/silicone	110 mm	98 mm	Barrel/plunger
5	10 ml PF JT	EVOC10PF	PP/PP/silicone	114 mm	104 mm	Barrel/plunger

ORAL SYRINGES ACCESSORIES



ORAL SYRINGE		1 ORAL SYRINGE WIPER (PP RESIN)	2 PROTECTIVE SHEATH (PP RESIN)	3 PIBA (PE RESIN) PP28 GLASS OR PLASTIC BOTTLE <i>CE marked</i>	4 CAP WITH PIBA INTEGRATED CRC CAP PP28 PIBA GLASS BOTTLE TE CAP DIN18 & CRC CAP DIN18 PIBA GLASS BOTTLE		5 SYRINGE HOLDER (PP RESIN)
Classic	0.3 ml						
	0.6 ml						X
	1.5 ml	X	X				X
	C2 ml	X					X
	2 ml		X				X
	TH5 ml						X
	5 ml	X	X				X
	C5 ml						X
	6 ml	X	X				X
	8 ml						
	10 ml						
Advanced	1 ml PF			X	X	X	
	3 ml PF			X	X	X	X
	TH5 ml LS						X
	TH5 ml PF			X	X	X	X
	C10 ml PF			X	X	X	
	C10 ml						
Evolution	1 ml PF JT			X	X	X	
	2 ml ORP						X
	3 ml PF JT			X	X	X	X
	TH5 ml PF JT			X	X	X	X
	10 ml PF JT			X	X	X	