

PRESS-IN-BOTTLE ADAPTOR PIBA



BENEFITS

- Flow restrictor that prevents accidental drug ingestion by children
- More hygienic: protects the medication from contamination during dispensing.
- No excess medication on the dosing syringe during dosing and administration, limiting the risk of overdose
- Allows the use of short oral dosing syringes with large bottles

PIBA	FOR BOTTLE USE	BOTTLE NECK	SYRINGE RANGES	SYRINGE TYPES
	PLASTIC	PP28	ADVANCED	FLAT END
	PLASTIC	PP28	ADVANCED / EVOLUTION	THIN TIP
	GLASS	PP28	ADVANCED	FLAT END
	GLASS	PP28	ADVANCED / EVOLUTION	THIN TIP

CHARACTERISTICS :

- CE marking
- Packaging and delivery options
 - Dosing syringe with reducer, assembled and non-wrapped
 - Dosing syringe with reducer, assembled and wrapped



ORAL SYRINGES ACCESSORIES



ORAL SYRINGE		1 ORAL SYRINGE WIPER (PP RESIN)	2 PROTECTIVE SHEATH (PP RESIN)	3 PIBA (PE RESIN) PP28 GLASS OR PLASTIC BOTTLE	4 CAP WITH PIBA INTEGRATED CRC CAP PP28 PIBA GLASS BOTTLE		5 SYRINGE HOLDER (PP RESIN)
Classic	0.3 ml						
	0.6 ml						X
	1.5 ml	X	X				X
	C2 ml	X					X
	2 ml		X				X
	TH5 ml						X
	5 ml	X	X				X
	C5 ml						X
	6 ml	X	X				X
	8 ml						
	10 ml						
Advanced	1 ml PF			X	X	X	
	3 ml PF			X	X	X	X
	TH5 ml LS						X
	TH5 ml PF			X	X	X	X
	C10 ml PF			X	X	X	
	C10 ml						
Evolution	1 ml PF JT			X	X	X	
	2 ml ORP						X
	3 ml PF JT			X	X	X	X
	TH5 ml PF JT			X	X	X	X
	10 ml PF JT			X	X	X	