





PRESS-IN-BOTTLE ADAPTOR **PIBA**



BENEFITS

- Flow restrictor that prevents accidental drug ingestion by children
- More hygienic: protects the medication from contamination during dispensing.
- No excess medication on the dosing syringe during dosing and administration, limiting the risk of overdose
- Allows the use of short oral dosing syringes with large bottles

PIBA	FOR BOTTLE USE	BOTTLE NECK	SYRINGE RANGES	SYRINGE TYPES
	PLASTIC	PP28	ADVANCED	FLAT END
	PLASTIC	PP28	ADVANCED / EVOLUTION	THIN TIP
	GLASS	PP28	ADVANCED	FLAT END
	GLASS	PP28	ADVANCED / EVOLUTION	THIN TIP

CHARACTERISTICS :

- CE marking
- Packaging and delivery options
 - Dosing syringe with reducer, assembled and non-wrapped
 - Dosing syringe with reducer, assembled and wrapped

